



Cambridge International AS & A Level

CANDIDATE
NAME

--

CENTRE
NUMBER

--	--	--	--	--	--

CANDIDATE
NUMBER

--	--	--	--

* 9 5 3 3 6 1 2 2 1 9 *



CHEMISTRY

9701/33

Paper 3 Advanced Practical Skills 1

May/June 2021

2 hours

You must answer on the question paper.

You will need: The materials and apparatus listed in the confidential instructions

INSTRUCTIONS

- Answer **all** questions.
- Use a black or dark blue pen. You may use an HB pencil for any diagrams or graphs.
- Write your name, centre number and candidate number in the boxes at the top of the page.
- Write your answer to each question in the space provided.
- Do **not** use an erasable pen or correction fluid.
- Do **not** write on any bar codes.
- You may use a calculator.
- You should show all your working, use appropriate units and use an appropriate number of significant figures.
- Give details of the practical session and laboratory, where appropriate, in the boxes provided.

Session
Laboratory

INFORMATION

- The total mark for this paper is 40.
- The number of marks for each question or part question is shown in brackets [].
- The Periodic Table is printed in the question paper.
- Notes for use in qualitative analysis are provided in the question paper.

For Examiner's Use	
1	
2	
3	
Total	

This document has **12** pages. Blank pages are indicated.

Quantitative analysis

Read through the whole method before starting any practical work. Where appropriate, prepare a table for your results in the space provided.


Show your working and appropriate significant figures in the final answer to **each** step of your calculations.

- 1 Bleach is made by reacting chlorine with a cold solution of sodium hydroxide. This reaction produces sodium chlorate(I), NaClO.




In this experiment you will determine the concentration of sodium chlorate(I) in a sample of bleach, **FA 1**.

To do this, you will react an acidified dilute solution of the bleach with iodide ions, I⁻. This reaction produces iodine, I₂.



The amount of iodine produced will then be determined by titration with thiosulfate ions, S₂O₃²⁻.



FA 1 is a solution of bleach.

FA 3 is dilute sulfuric acid, H₂SO₄.

FA 4 is 0.500 mol dm⁻³ potassium iodide, KI.

FA 5 is 0.100 mol dm⁻³ sodium thiosulfate, Na₂S₂O₃.
starch indicator

(a) Method

Dilution

- Pipette 25.0 cm³ of **FA 1** into the 250 cm³ volumetric flask.
- Add distilled water to make 250 cm³ of solution and shake the flask thoroughly.
- Label this flask **FA 2**.

Titration

- Fill a burette with **FA 5**.
- Rinse the pipette thoroughly with distilled water and then with a little **FA 2**.
- Pipette 25.0 cm³ of **FA 2** into a conical flask.
- Use the measuring cylinder to add 20 cm³ of **FA 3** to the conical flask.
- Use the measuring cylinder to add 15 cm³ of **FA 4** to the conical flask. The solution will turn brown as iodine is produced.
- Add **FA 5** from the burette until the solution has turned yellow.
- Add 10 drops of starch indicator to the conical flask. The solution will turn blue-black.
- Continue to add more **FA 5** from the burette until the blue-black colour just disappears. This is the end-point of the titration.
- Carry out a rough titration and record your burette readings in the space provided.

The rough titre = cm³.

- Carry out as many accurate titrations as you think necessary to obtain consistent results.
- Make sure that your recorded results show the precision of your practical work.
- Record in a suitable form in the space below all your burette readings and the volume of **FA 5** added in each accurate titration.

Keep FA 4 and FA 5 for use in Question 3.

I	
II	
III	
IV	
V	
VI	
VII	

[7]

- (b)** From your accurate titration results, obtain a value for the volume of **FA 5** to be used in your calculations. Show clearly how you obtained this value.

25.0 cm³ of **FA 2** required cm³ of **FA 5**.
[1]

(c) Calculations

(i) Give your answers to **(c)(ii)** and **(c)(iii)** to the appropriate number of significant figures. [1]

(ii) Use your answer to **(b)** and the relevant equation on page 2 to calculate the number of moles of iodine that formed when 25.0 cm³ of **FA 2** reacted with **FA 4**.

moles of I₂ = mol [1]

(iii) Calculate the concentration of sodium chlorate(I) in **FA 1**. Show your working.

concentration of NaClO = mol dm⁻³ [2]

(d) (i) In this method an excess of potassium iodide must be added to 25.0 cm³ of **FA 2**.

Use your answer to **(c)(ii)** to show by calculation that the potassium iodide is in excess.

[1]

(ii) A student carries out the same method but the concentration of the potassium iodide solution is not stated.

What change to the practical procedure could the student make to check that the potassium iodide is in excess?
Explain your answer.

.....
.....
..... [2]

[Total: 15]

- 2 In this experiment you will determine the value of x in the formula for hydrated sodium thiosulfate, $\text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O}$. In this formula x is an integer.

FA 6 is hydrated sodium thiosulfate, $\text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O}$.

(a) Method

- Support the cup in the 250 cm^3 beaker.
- Rinse the measuring cylinder thoroughly with distilled water.
- Use the measuring cylinder to transfer 20.0 cm^3 of distilled water into the cup.
- Measure and record the initial temperature of the water in the cup. Tilt the cup if necessary so that the bulb of the thermometer is fully covered.
- Weigh the stoppered container of **FA 6** and record the mass.
- Carefully add all the sample of **FA 6** to the water in the cup.
- Stir the mixture and record the minimum temperature.
- Reweigh the stoppered container and any residual **FA 6**. Record the mass.
- Calculate and record the mass of **FA 6** added to the water and the change in temperature.

I	
II	
III	
IV	

[4]

(b) Calculations

- (i) Calculate the heat energy change in the reaction.
(Assume that 4.2 J are required to change the temperature of 1.0 cm^3 of solution by $1.0\text{ }^\circ\text{C}$.)

heat energy change = J
[1]

- (ii) Use your answer to (b)(i) to calculate the relative formula mass, M_r , of hydrated sodium thiosulfate.

You should assume that the enthalpy change of hydration, ΔH_{hyd} , for $\text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O}$ is $+47.4 \text{ kJ mol}^{-1}$ under the conditions of your experiment.

Show your working.

$$M_r \text{ of } \text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O} = \dots\dots\dots [2]$$

- (iii) Use your answer to (b)(ii) to calculate the value of x in the formula $\text{Na}_2\text{S}_2\text{O}_3 \cdot x\text{H}_2\text{O}$.

$$x = \dots\dots\dots [1]$$

- (c) A student repeats the experiment but adds 25.0 cm^3 of distilled water into the cup instead of 20.0 cm^3 . The student carries out the calculations based on adding 20.0 cm^3 .

What effect would this have on the student's answer to (b)(ii)?

Explain your answer.

.....

.....

..... [2]

[Total: 10]

Qualitative analysis

Where reagents are selected for use in a test, the **name** or **correct formula** of the element or compound must be given.

At each stage of any test you are to record details of the following:

- colour changes seen
- the formation of any precipitate and its solubility in an excess of the reagent added
- the formation of any gas and its identification by a suitable test.

You should indicate clearly at what stage in a test a change occurs.

If any solution is warmed, a **boiling tube** must be used.

Rinse and reuse test-tubes and boiling tubes where possible.

No additional tests for ions present should be attempted.

3 (a) FA 7 is a salt containing one cation and one anion, both of which are listed in the Qualitative Analysis Notes.

- (i) Transfer a small spatula measure of **FA 7** to a hard-glass test-tube and heat, gently at first and then more strongly.
Record all your observations.

.....
.....
..... [2]

- (ii) Place the remaining **FA 7** in the 100 cm³ beaker and dissolve in approximately 50 cm³ of distilled water.
Use a sample of this solution to test for the cation present in **FA 7**.
Record your tests and observations in the space below.

Keep the solution for use in (a)(iii).

[3]

- (iii) To a 1 cm depth of solution **FA 7** in a test-tube, add a few drops of aqueous silver nitrate. Record your observations.

..... [1]

- (iv) Suggest the ions present in **FA 7**.

cation anion [1]

- (b) **FA 8** is acidified aqueous iron(III) chloride, FeCl_3 .

Carry out the following tests and record your observations.

<i>test</i>	<i>observations</i>
<p>Test 1 To a 1 cm depth of FA 8 in a test-tube, add a 1 cm depth of FA 4, then</p>	
<p>-----</p> <p>add starch indicator.</p>	
<p>Test 2 To a 1 cm depth of FA 8 in a test-tube, add a 1 cm depth of FA 5. Allow to stand until there is no further change, then</p>	
<p>-----</p> <p>add a 1 cm depth of FA 4, then</p>	
<p>-----</p> <p>add starch indicator.</p>	

[3]

- (c) Half fill the 250 cm³ beaker with water. Place the beaker on a gauze supported on the tripod and heat the water to approximately 70 °C. Switch off the Bunsen burner. This will be used as a water-bath.

FA 9, FA 10 and FA 11 are each one of the following.

- 1.0 mol dm⁻³ hydrochloric acid, HCl
- 1.0 mol dm⁻³ methanoic acid, HCOOH
- 1.0 mol dm⁻³ sodium sulfite, Na₂SO₃

- (i) Carry out the following tests and record your observations. For each test use a 1 cm depth in a test-tube.

test	observations		
	FA 9	FA 10	FA 11
Test 1 Add a 1 cm strip of magnesium.			
Test 2 Add a few drops of aqueous acidified potassium manganate(VII). Place the test-tube in the hot water-bath.			

[3]

- (ii) Identify each of the solutions.

FA 9 is

FA 10 is

FA 11 is

[2]

[Total: 15]

Qualitative analysis notes

1 Reactions of aqueous cations

ion	reaction with	
	NaOH(aq)	NH ₃ (aq)
aluminium, Al ³⁺ (aq)	white ppt. soluble in excess	white ppt. insoluble in excess
ammonium, NH ₄ ⁺ (aq)	no ppt. ammonia produced on heating	–
barium, Ba ²⁺ (aq)	faint white ppt. is nearly always observed unless reagents are pure	no ppt.
calcium, Ca ²⁺ (aq)	white ppt. with high [Ca ²⁺ (aq)]	no ppt.
chromium(III), Cr ³⁺ (aq)	grey-green ppt. soluble in excess	grey-green ppt. insoluble in excess
copper(II), Cu ²⁺ (aq)	pale blue ppt. insoluble in excess	pale blue ppt. soluble in excess giving dark blue solution
iron(II), Fe ²⁺ (aq)	green ppt. turning brown on contact with air insoluble in excess	green ppt. turning brown on contact with air insoluble in excess
iron(III), Fe ³⁺ (aq)	red-brown ppt. insoluble in excess	red-brown ppt. insoluble in excess
magnesium, Mg ²⁺ (aq)	white ppt. insoluble in excess	white ppt. insoluble in excess
manganese(II), Mn ²⁺ (aq)	off-white ppt. rapidly turning brown on contact with air insoluble in excess	off-white ppt. rapidly turning brown on contact with air insoluble in excess
zinc, Zn ²⁺ (aq)	white ppt. soluble in excess	white ppt. soluble in excess

2 Reactions of anions

<i>ion</i>	<i>reaction</i>
carbonate, CO_3^{2-}	CO_2 liberated by dilute acids
chloride, $\text{Cl}^-(\text{aq})$	gives white ppt. with $\text{Ag}^+(\text{aq})$ (soluble in $\text{NH}_3(\text{aq})$)
bromide, $\text{Br}^-(\text{aq})$	gives cream ppt. with $\text{Ag}^+(\text{aq})$ (partially soluble in $\text{NH}_3(\text{aq})$)
iodide, $\text{I}^-(\text{aq})$	gives yellow ppt. with $\text{Ag}^+(\text{aq})$ (insoluble in $\text{NH}_3(\text{aq})$)
nitrate, $\text{NO}_3^-(\text{aq})$	NH_3 liberated on heating with $\text{OH}^-(\text{aq})$ and Al foil
nitrite, $\text{NO}_2^-(\text{aq})$	NH_3 liberated on heating with $\text{OH}^-(\text{aq})$ and Al foil
sulfate, $\text{SO}_4^{2-}(\text{aq})$	gives white ppt. with $\text{Ba}^{2+}(\text{aq})$ (insoluble in excess dilute strong acids)
sulfite, $\text{SO}_3^{2-}(\text{aq})$	gives white ppt. with $\text{Ba}^{2+}(\text{aq})$ (soluble in excess dilute strong acids)

3 Tests for gases

<i>gas</i>	<i>test and test result</i>
ammonia, NH_3	turns damp red litmus paper blue
carbon dioxide, CO_2	gives a white ppt. with limewater (ppt. dissolves with excess CO_2)
chlorine, Cl_2	bleaches damp litmus paper
hydrogen, H_2	'pops' with a lighted splint
oxygen, O_2	relights a glowing splint

The Periodic Table of Elements

		Group																																	
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18																		
		<table border="1" style="margin: auto;"> <tr> <td style="text-align: center;">1</td> <td style="text-align: center;">H hydrogen 1.0</td> </tr> </table>																1	H hydrogen 1.0																
1	H hydrogen 1.0																																		
		<table border="1" style="margin: auto;"> <tr> <td style="text-align: center;">atomic number</td> <td style="text-align: center;">atomic symbol</td> <td style="text-align: center;">name</td> <td style="text-align: center;">relative atomic mass</td> </tr> </table>																atomic number	atomic symbol	name	relative atomic mass														
atomic number	atomic symbol	name	relative atomic mass																																
3	4	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32	33	34	35	36								
Li lithium 6.9	Be beryllium 9.0	Na sodium 23.0	Mg magnesium 24.3	Al aluminium 27.0	Si silicon 28.1	P phosphorus 31.0	S sulfur 32.1	Cl chlorine 35.5	Ar argon 39.9	K potassium 39.1	Ca calcium 40.1	Sc scandium 45.0	Ti titanium 47.9	V vanadium 50.9	Cr chromium 52.0	Mn manganese 54.9	Fe iron 55.8	Co cobalt 58.9	Ni nickel 58.7	Cu copper 63.5	Zn zinc 65.4	Ga gallium 69.7	Ge germanium 72.6	As arsenic 74.9	Se selenium 79.0	Br bromine 79.9	Kr krypton 83.8								
37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57-71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86
Rb rubidium 85.5	Sr strontium 87.6	Y yttrium 88.9	Zr zirconium 91.2	Nb niobium 92.9	Mo molybdenum 95.9	Tc technetium —	Ru ruthenium 101.1	Rh rhodium 102.9	Pd palladium 106.4	Ag silver 107.9	Cd cadmium 112.4	In indium 114.8	Sn tin 118.7	Pb lead 207.2	Bi bismuth 209.0	Po polonium —	At astatine —	Xe xenon 131.3	Cs caesium 132.9	Ba barium 137.3	lanthanoids	Hf hafnium 178.5	Ta tantalum 180.9	W tungsten 183.8	Re rhenium 186.2	Os osmium 190.2	Ir iridium 192.2	Au gold 197.0	Hg mercury 200.6	Tl thallium 204.4	Pb lead 207.2	Bi bismuth 209.0	Po polonium —	At astatine —	Rn radon —
87	88	89-103	104	105	106	107	108	109	110	111	112	113	114	115	116	117	118	119	120	121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136
Fr francium —	Ra radium —	actinoids	Rf rutherfordium —	Db dubnium —	Sg seaborgium —	Bh bohrium —	Hs hassium —	Mt meitnerium —	Ds darmstadtium —	Rg roentgenium —	Cn copernicium —	Nh nihonium —	Fl flerovium —	Lv livermorium —	Ts tennessine —	Og oganesson —	118	119	120	121	122	123	124	125	126	127	128	129	130	131	132	133	134	135	136

lanthanoids	57	58	59	60	61	62	63	64	65	66	67	68	69	70	71
	La lanthanum 138.9	Ce cerium 140.1	Pr praseodymium 140.9	Nd neodymium 144.4	Pm promethium —	Sm samarium 150.4	Eu europium 152.0	Gd gadolinium 157.3	Tb terbium 158.9	Dy dysprosium 162.5	Ho holmium 164.9	Er erbium 167.3	Tm thulium 168.9	Yb ytterbium 173.1	Lu lutetium 175.0
actinoids	89	90	91	92	93	94	95	96	97	98	99	100	101	102	103
	Ac actinium —	Th thorium 232.0	Pa protactinium 231.0	U uranium 238.0	Np neptunium —	Pu plutonium —	Am americium —	Cm curium —	Bk berkelium —	Cf californium —	Es einsteinium —	Fm fermium —	Md mendelevium —	No nobelium —	Lr lawrencium —